



Diffusion with Pulsed NMR


Kyle Moe

Introduction




- Nuclear Magnetic Resonance (NMR)
 - Place atomic nuclei into background magnetic field B_0
 - Apply circularly polarized perturbation field B_1
 - On resonance, nuclei precess around B_1
- Applications of this phenomenon
 - Magnetic Resonance Imaging (MRI)
 - Material properties
 - Rate of alignment with background magnetic field (T_1)
 - Rate of magnetization decay (T_2)
 - Diffusion coefficient (D) of material
 - Spectroscopy to reveal molecular structure

Protons in B_0



$$B_0 = 0$$

$$B_0 \neq 0$$



$$\mu_z = \gamma m_I = \pm \gamma \hbar / 2$$

$$U = -\vec{\mu} \cdot \vec{B}$$

$$\Delta U = \gamma \hbar B_0$$


$$\Delta U = \hbar \omega_0$$

$$\omega_0 = \gamma B_0$$


NMR Apply Perturbation Field B_1

$$\mathbf{B}_{\text{eff}}^* = B_1 \hat{\mathbf{i}}^* + (B_0 - \omega/\gamma) \hat{\mathbf{k}}^*$$

- Perpendicular to B_0
- Circularly polarized
- Oscillating at ω




In rotating frame



$\pi/2$ pulse $\rightarrow M_z \rightarrow M_{x,y}$
 π pulse $\rightarrow M_z \rightarrow -M_z$
 2π pulse $\rightarrow M_z \rightarrow M_z$

In stationary frame




Measuring T_2 and D

Measuring the decay of spin echoes


Employ the following pulse sequence:

- Flip the spins into the x-y plane with $\pi/2$ pulse
- Retain spins in x-y plane with n number of π pulses



$$M_{x,y}(t) = M_0 e^{-\left(\frac{1}{T_2} - \frac{\gamma^2 \left(\frac{\partial B_{0,z}}{\partial x,y}\right)^2 D t^2}{12n^2}\right)t}$$

In rotating frame






Limitations of Experiment/Apparatus



- Highly non-uniform magnetic field
 - Accurate T_2 measurements require suppression of Diffusion term
 - Gradient changes across width of sample
- Apparatus constraints
 - Maximum number of π pulses (99)
 - Limited frequency range of detector coil


Experimental Apparatus



Exploring Permanent Magnet

On resonance: $\omega = \gamma B_0$

- Tune TeachSpin NMR apparatus to this frequency
 - Use this frequency to derive B_0
- What is the utility of this?
 - Identify low-gradient regions
 - Identify regions of near-linear gradient



What does tuning to resonance mean?

Aliased frequency between larmor and applied frequencies



vs.

Aliased frequency approaches 0Hz on resonance



Measuring Decay (T_2)

Measurements made of echo peaks (V)



Delay time between peaks gives time (t)

Using a π -pulse sequence with many pulses, we can measure T_2 by assuming:

$$M_{x,y}(t) = M_0 e^{-t/T_2}$$

- Requires that Diffusion is suppressed
 - Low gradient region
 - High density of pulses


~~$$M_{x,y}(t) = M_0 e^{-\left(\frac{1}{T_2} - \frac{\gamma^2 \left(\frac{\partial B_{0,z}}{\partial x,y}\right)^2 D t^2}{12 n^2}\right) t}$$~~



Measuring Diffusion Coefficient (D)

- Similar pulse sequence
 - Fewer π -pulses
 - Larger gradient
- With known T_2 we can fit with D as the only unknown parameter

$$M_{x,y}(t) = M_0 e^{-\left(\frac{1}{T_2} - \frac{\gamma^2 \left(\frac{\partial B_{0,z}}{\partial x,y}\right)^2 D t^2}{12 n^2}\right) t}$$




Magnetic Field Gradient Measurement

Assume T_2 measurement region has gradient of 0.

Gradient in diffusion region:


$$\left| \frac{\partial B_{0,z}}{\partial x, y} \right| = 2.6 \pm 0.2 \text{ gauss/cm}$$

This value is used for all future diffusion measurements




T₂ and D of Water


Water has both a long T₂ and high D



vs.



Decay of Magnetization of H₂O with and without Diffusion



$$M_{x,y}(t) = M_0 e^{-\left(\frac{1}{T_2} - \frac{\gamma^2 \left(\frac{\partial B_{0,z}}{\partial x,y}\right)^2 D t^2}{12 n^2}\right) t}$$

$$T_2 = 2.58 \pm 0.03 \text{ s}$$

$$D = (2.0 \pm 0.3) * 10^{-5} \text{ cm}^2 \text{ s}^{-1}$$

(at 21°C)

Most of the uncertainty in D comes from uncertainty in the gradient

Substances of Complex Molecules

- Much thicker than water or isopropanol.
 - High viscosity means it is harder for molecules to diffuse through to regions of different magnetic field
- Decays are more complex
 - Biexponential*

$$M_{x,y}(t) = M_0 e^{-t/T_2}$$




$$M_{x,y}(t) = M_{0,1} e^{-t/T_{2,1}} + M_{0,2} e^{-t/T_{2,2}}$$

$$\text{where } M_{0,1} + M_{0,2} = M_0$$

*May be additional terms that are less prevalent




3 in 1 Household Oil



Light Mineral Oil

Mineral Oil T_2 and D

$$M_{x,y}(t) = (M_{0,1}e^{-\frac{t}{T_{2,1}}} + M_{0,2}e^{-\frac{t}{T_{2,2}}})e^{-\frac{\gamma^2(\frac{\partial B_{0,z}}{\partial x,y})^2 Dt^3}{12n^2}}$$



$$T_{2,1} = 9.4 \pm 0.1 \text{ ms}$$


$$T_{2,2} = 49.1 \pm 0.5 \text{ ms}$$

$$D = (2.5 \pm 0.6) * 10^{-7} \text{ cm}^2 \text{ s}^{-1}$$

D is much smaller than that of water.

3 in 1 T_2 and D

$$M_{x,y}(t) = (M_{0,1}e^{-\frac{t}{T_{2,1}}} + M_{0,2}e^{-\frac{t}{T_{2,2}}})e^{-\frac{\gamma^2(\frac{\partial B_{0,z}}{\partial x,y})^2 Dt^3}{12n^2}}$$



$$T_{2,1} = 29 \pm 1 \text{ ms}$$


$$T_{2,2} = 131 \pm 5 \text{ ms}$$

$$D = (4.8 \pm 0.8) * 10^{-7} \text{ cm}^2 \text{ s}^{-1}$$

Effect of Diffusion




Decays with and without Diffusion 3 in 1



$$D = (4.8 \pm 0.8) * 10^{-7} \text{ cm}^2 \text{ s}^{-1}$$

Decay of Magnetization of H2O with and without Diffusion




$$D = (2.0 \pm 0.3) * 10^{-5} \text{ cm}^2 \text{ s}^{-1}$$

NMR Spectroscopy on Mineral Oil




Free-Induction Decay





Discrete
Fourier
Transform



Full-Range Discrete Fourier Transform of FID



Existence of Distinct Peaks



NMR as a Method of Measuring Diffusion



- From the results of this experiment, we can broadly say that NMR provides an effective means of measuring D
 - D of water was consistent
 - Low D -> less stark effect of diffusion
- Not all protons are the same
 - Need for bi-exponential fit for mineral oil and household oil
 - Local environments change interaction
 - Smaller molecules have longer T_2 from molecular tumbling

Further Work



- Identify exactly which protons in mineral oil contribute to each decay
- Explore broader range of materials
- Refine measurements of T_2 by reducing size of sample and find region of near-zero gradient
- Reduced sample size allows for linear gradient approximation to be more accurate as well
- Isopropanol is complex and may be characterized by more than one T_2
- D changes with temperature